

# EarGenie - A clinician-friendly fNIRs system to evaluate infant speech detection & discrimination

G. Balasubramanian<sup>1,2</sup>, J. Wunderlich<sup>1,2</sup>, D. Mao<sup>1,2</sup>, L. McDonald<sup>1</sup>, D. X. Gao<sup>1,2</sup>, C. M. McKay<sup>1,2</sup>

<sup>1</sup> The Bionics Institute, 384-388 Albert St, East Melbourne VIC 3002, Australia, <sup>2</sup> Department of Medical Bionics, University of Melbourne, Parkville, VIC 3010, Australia.

\* For more information, please email [GBalasubramanian@bionicsinstitute.org](mailto:GBalasubramanian@bionicsinstitute.org)

## Introduction

### Infant Speech Detection & Discrimination

- Current methods of auditory testing cannot accurately test for speech discrimination which is crucial for language development.
- Hearing assessments for infants with auditory neuropathy (AN) are also not reliable.
- Our system, EarGenie (Fig. 1), is an end-to-end fNIRs (Functional Near Infra Red Spectroscopy) system that tests for detection and discrimination including for infants with AN.
- It can be used to evaluate hearing aid and Cochlear Implant (CI) programs due to its patient-friendly, non-invasive nature.
- Our custom analysis algorithm was developed and evaluated using NIRx sensor measurements across 4 regions of interest (ROI): left and right temporal and prefrontal regions.

### EarGenie - Prototype build in progress




Fig. 1: Ear Genie System components

## Method

### Participants & Stimuli

- We show results for speech detection ("Ba" presented at 65 dB SPL) and discrimination using 3 different speech-token contrasts: Ba/Tea, Ba/Bee, Ba/Ga as well as performance across tokens presented at multiple intensities.
- Subjects were 34 infants with normal hearing aged between 2 and 20 months who were tested while asleep. The stimulus presentation protocol used is shown in Fig 2.
- Each 5.4 s stimuli block, was either 12 concatenated "Ba" tokens in the detection phase or alternating "Ba" and contrast token in the discrimination phase.




Fig. 2: Stimulus presentation paradigm. To measure speech detection, the standard blocks were separated by silence; to measure speech discrimination, the novel stimuli blocks were separated by a repeated Ba stimulus.

### Data Pre-processing included:

- Conversion to optic density and removal of noisy channels and physiological noise using the Scalp Coupling Index (>0.75), Temporal Derivative Dispersion Repair and Bandpass Filter (0.01 – 0.25 Hz).
- Conversion to HbO and HbR signals by applying modified Beer Lambert law.
- Epoching and baseline correction of fNIRs signals between -3 and 27 s post stimulus onset. Control epochs were also gathered from the silence (detection) or non-silence baseline (discrimination).

### Post-Processing:

- The EarGenie custom algorithm uses online processing. Hence data was processed in an incremental fashion, starting after 3 trials, and after each subsequent trial. The testing automatically stops the test for a token when it identifies a significant detection or discrimination.
- A stochastic process approximates epochs gathered in both control and post-stimulus conditions.
- Statistical tests are carried out to evaluate the similarity of control and stimulus epochs. The statistical tests are permuted by varying the controls to ensure the presence of a consistent response.
- A control vs. control test is executed to test for the possibility of a false positive. Additionally, separate testing was conducted with silent 0 dB trigger in a subset of subjects.

## Results

We require our algorithm to show a significant result in at least 2 of the 4 ROIs to be considered a true detection or discrimination. Tables 1 and 2 show the number of infants out of 32 who showed significant results in different number of ROIs.

Table 1 - Speech Detection Sensitivity

# ROIs Detected out of 4	All 4 ROIs	3 ROIs	2 ROIs	1 ROIs	0 ROIs
# infants out of 32	21	9	1	1	0

- For detection (Table 1), 31 infants out of 32 showed detection in 2 or more ROIs (97 % sensitivity).

Table 2 - Speech Discrimination Sensitivity

# ROIs Detected out of 4	All 4 ROIs	3 ROIs	2 ROIs	1 ROIs	0 ROIs
Ba/Tea # infants out of 16	8	5	2	1	0
Ba/Bee # infants out of 26	15	5	3	1	2
Ba/Ga # infants out of 18	9	3	2	3	1

- For discrimination (Table 2), 15 infants of 16 showed Ba/Tea discrimination (94 % sensitivity).
- 23 infants out of 26 showed Ba/Bee discrimination (88 % sensitivity).
- 14 infants out of 18 showed Ba/Ga discrimination (78 % sensitivity).

Table 3 - Speech Detection Sensitivity at soft intensity levels

# ROIs Detected out of 4	≥ 2 ROIs	< 2 ROIs
35 dB SPL # infants out of 6	6	0
45 – 50 dB SPL # infants out of 13	11	2

- For detection at soft intensity levels (Table 3), 6 infants out of 6 showed detection at 35 dB SPL in 2 or more ROIs (100 % sensitivity).
- 11 infants out of 13 showed detection at 45 - 50 dB SPL in 2 or more ROIs (85 % sensitivity).

### Specificity

The control vs. control tests were conducted across all subjects and speech tokens to identify the probability of false positives. The probability of false positives is 5 %, i.e., 95% specificity.

### Test efficiency

A significant detection or discrimination response was identified by the algorithm within 5 trials (equivalent to 2.5 minutes following the control segment) in ~ 70 % of the cases described, which validates an online processing approach.

## Conclusion

We have a robust algorithm to determine detection and discrimination of speech sounds by sleeping infants using fNIRs sensors. Furthermore, significant response detection was often obtained using less than 5 sound presentations, and at levels as low as 35 dB SPL. Further testing with the EarGenie hardware is set to commence shortly after bench testing success.

### References:

- [1] McKay et. al. A reliable, accurate, and clinician-friendly objective test of speech sound detection and discrimination in sleeping infants. (PsyArxiv)
- [2] McDonald et. al. EarGenie, an innovative test to measure speech discrimination using functional near-infrared spectroscopy. J. Acoust. Soc. Am. 154, A33-A34 (2023)

## Acknowledgements

The authors would like to thank the Victorian Infant Hearing Screening Program (VIHSP) for its support in participant recruitment. This project was supported by The National Health and Medical Research Council of Australia, the Australian Government through MTPConnect and the Medical Research Future Fund, the Victoria State Government through the Victorian Medical Research Acceleration Fund, and the Garnett Passe and Rodney Williams Memorial Foundation. The Bionics Institute acknowledges the support it receives from the Victorian State Government through its Operational Infrastructure Support Program.